**AGC®**  
Galvanotechnik

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AGC®

electroforming technique

AGC® Manual  
for Dentists and Dental Technicians

english

■-0483

**WIELAND**  
DENTAL • DENTURALS

■ Thank you for choosing the  
AGC® electroforming technique from  
Wieland Dental + Technik.

The AGC® process has been successfully implemented for dental restorations since 1986.

Today it is one of the most tested and future-oriented technologies for producing high-quality crowns and bridges. The AGC® electroforming technique has been constantly subjected to scientific studies carried out by reputable universities. Take a little time to read through these instructions carefully before you first use the equipment. This will help to avoid incorrect operation, which could lead to unsatisfactory results.

Should you have any questions after studying these instructions, please contact your distributor or directly:

Wieland Dental + Technik GmbH & Co. KG  
Schwenninger Straße 13  
75179 Pforzheim/Germany  
Tel.: +49 72 31/37 05-0  
Fax: +49 72 31/35 79 59  
Internet: <http://www.wieland-dental.de>  
eMail: [info@wieland-dental.de](mailto:info@wieland-dental.de)

■ Technical data:  
AGC® electroformed gold

|                                                   |                                                                                                                |
|---------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| Composition                                       | > 99.9 %                                                                                                       |
| Color                                             | golden-yellow                                                                                                  |
| Biologically tested                               | ✓                                                                                                              |
| Density, g/cm <sup>3</sup>                        | 19.3                                                                                                           |
| Melting point °C                                  | 1064                                                                                                           |
| Vickers hardness,<br>HV 0.1/30                    | n 30 g > 120                                                                                                   |
| CTE (25-500 °C), 10 <sup>-6</sup> K <sup>-1</sup> | 15.5                                                                                                           |
| Elasticity modulus, MPa                           | 80.000                                                                                                         |
| Suitable<br>dental ceramic                        | Firing temperature<br>max. 950 °C<br>CTE (25-500 °C)<br>approx. 11 to<br>16 · 10 <sup>-6</sup> K <sup>-1</sup> |

n = after firing

g = after electroforming

# ■ Contents

Areas of Application for  
AGC® Electroforming Units ..... p. 4

Dental procedures

General introduction ..... p. 5

Dental preparations ..... p. 6

Taking the impression ..... p. 7

Fitting and cementation ..... p. 8

Various uses and fitting  
of AGC® SupraCaps® ..... p. 9

Technical procedures

Crowns and anchor crowns ..... p. 10

§Model making/duplicating ..... p. 10

§Duplicating the die ..... p. 11

§Microwave drying ..... p. 12

§Contact rods/conductive silver lacquer ..... p. 12

§Margin correction ..... p. 14

§Veneering with ceramic and resin ..... p. 14

AGC® electroformed bridge –  
sintering procedure ..... p. 16

AGC® SupraCaps® on implant posts ..... p. 19

AGC® telescope crown technique ..... p. 21

AGC® Bar sleeves ..... p. 26

## ■ Areas of Application for AGC® Electroforming Units

| Indications                                   | AGC® 5 Process | AGC® Micro | AGC® Micro Plus | AGC® Speed |
|-----------------------------------------------|----------------|------------|-----------------|------------|
| Crowns (A) (veneered)                         | ✓              | ✓          | ✓               | ✓          |
| Splints (A)                                   | ✓              | ✓          | ✓               | (✓)*       |
| Anchor crowns<br>for bridges (C)              | ✓              | ✓          | ✓               | ✓          |
| Electroformed<br>small bridges (C)            | ✓              | ✓          | ✓               |            |
| Secondary crowns<br>(Telescope technique) (B) | ✓              | ✓          | ✓               | ✓          |
| SupraCaps®<br>Implant prosthetics (D)         | ✓              | ✓          | ✓               | ✓          |
|                                               | Order-No.      | Order-No.  | Order-No.       | Order-No.  |
| Use AGC®<br>gold electrolyte (CE-0483)        | 6622           | 6607       | 6607            | 6850       |
| Use AGC®<br>Brightener (CE-0483)              | 6674           | 6674       | 6674            | 6851       |

\* max. 1.15 g

AGC® electroformed substructures must always be built-up with ceramic or composites.

# Dental procedures

## ■ General introduction

When supplying patients with electroformed restorations, the usual rules of dental prosthetics still apply.

Patient selection and preparation should take place *lege artis*.

The selection and prognosis of suitable abutment teeth is as for conventional procedures.

## ■ Indications

The electroforming technique can be used for the following prosthetic restorations:

- A Veneered single crowns
- B Combination of fixed and removable restorations: telescope crown technique
- C Small bridges: One pontic in the upper anterior region and upper and lower posterior region; two pontics in the lower anterior region.
- D Implant suprastructures

## ■ Contraindications

Bruxism or traumatic occlusal relations.

Too little space for electroformed crowns built-up with ceramic, eg. antagonist is too close. (The tooth preparation must be anatomical with a reduction of 1.2 – 1.4 mm.)

Clinically old crowns which no longer allow for *lege artis* preparation.

Tangential preparation (the preferred type of preparation is a circumferential chamfer or shoulder with rounded inner edges).

Aptyalism (saliva deficiency) can cause functional problems in the telescope crown technique.

Telescope primary crowns containing copper can cause the inner surfaces of the AGC® secondary crown to discolor.

Temporary insertion of electroformed restorations in contra-indicated.


Specific technical requirements and special characteristics, with regard to the different indication classes, are referred to in the technical procedures. (Restrictions regarding the manufacture of the indication classes may apply in the case of certain types of units; see the chart entitled "Areas of Application for AGC® Electroforming Units".)

The following description refers in general to all indication classes and deals exclusively with the special characteristics and requirements applicable to the preparation, implementation and fitting of dental restorations fabricated by the AGC® electroforming technique.

## ■ Dental preparations

### ■ Preparation of abutment teeth

The abutment tooth must be built up and filled prior to the preparation process. Subsequent blocking out performed on the model or by using cementing material should be avoided, firstly because any areas of abutment tooth which are not blocked out can lead to differences of dimension and changes of shape during the subsequent production of the model and secondly because cemented grooves which are too thick can lead to a loosening of the cement layer and cause the restoration to come away.

### ■ Crowns and anchor crowns

Any tooth which requires covering with a crown can benefit from a long-term restoration in the form of an electroformed crown. This statement is made on the assumption that the practitioner exercises the required care both during the preparation and impression-taking as during fitting of the restoration. In order to receive an electroformed crown the tooth die must be prepared in such a way that it is as smooth as possible, without ridges or other irregularities and with a clear edge to the preparation. The preferred type of preparation is a clear chamfer; a **tangential preparation is contra-indicated**.

The tooth can be prepared using a torpedo bur ISO size 014. The pre-preparation can be carried out using coarse diamond burs (green ring). The edge of the preparation should end at least 1 mm above the gingival margin. The edge of the preparation can be lowered in two stages, either first with a torpedo bur of the same size and with a red ring (fine-grain) or with a torpedo bur of the same type, but with a yellow ring (15 µm grain). The preparation must follow the anatomic form of the abutment tooth.

The aim of the preparation is to reduce the tooth substance both around its circumference and in the occlusal plane by 1.2 – 1.4 mm, in order to provide enough room for a **ceramic layer of 1 – 1.2 mm**.

This is especially important for the incisal edge and the molar cusps. If too much of the cusp is taken away ( $\geq 1.7\text{mm}$ ), then too much ceramic will be built up which subsequently causes a higher risk of fracture.

Optionally, cylindrical preparation diamonds (stepped preparation with rounded cylindrical bur ISO size 0.12/0.14) can be used. In this case it is important to avoid the danger of loss of vitality resulting from overpreparation (particularly in the case of mandibular front teeth, maxillary lateral incisors and maxillary pre-molars).

**Undercutting areas of preparation which have to be blocked out in the laboratory are contra-indicated and must similarly be avoided.**

## ■ Taking the impression

Because of the accuracy and precision of the AGC® electroforming process, the impression must be taken with particular care. This applies equally to crowns as well as anchor crowns. It is important for the impression to reproduce the edge of the preparation in an area apical to the preparation edge exactly, so that it can subsequently be transferred to the situation on the model.

### ■ Gingival retraction

In the case of epigingival and subgingival preparation edges it is particularly difficult to transfer the situation in the mouth to that on the model. If the patient shows a good state of parodontal health, then the edge of the preparation can be reproduced accurately lege artis by using the method of gingival retraction. For this purpose retraction threads or rings can be used. Depending on the depth of sulcus and the clinical condition these threads and rings can be combined in a variety of thickness and layers. The process of gingival retraction must

be carried out with great care. Retraction threads or rings must be laid carefully and must not be inserted with force. A traumatic procedure results in the destruction of the fibre structure causing subsequent inflammation and osseoresorption leading to unnecessary recession once the margin has been restored to its biological width.

### ■ Impression technique

Impressions can be taken with all commercially available elastomer impression materials. Reversible elastic hydrocolloid impression materials may also be used. When using elastomers it is advisable to use customised impression trays. The recommendations of the manufacturer of the impression material must be followed at all times. It is easier to reproduce the details when preparations have a supra-gingival taper. However, both clinical situations require a standardised procedure for taking the impression. Firstly the abutment teeth must be relatively dry for the impression to be taken and secondly, care must be taken to remove blood, saliva and tooth grindings in preparation. It is impossible to avoid squirting impression material onto the abutment tooth. It may even be necessary to blast away material from the die with the aid of an air jet. When it is applied, the impression material must be free of cavitations both around the abutment tooth and in the impression tray. After setting, the impression tray is removed from the mouth and the impression is checked for defects. It is important to ensure that no part of the impression is penetrated in the tooth area. The edges of the preparation at the abutment tooth must be fully reproduced in the impression. The model is then made in the laboratory.

## ■ Fitting and cementation

### ■ Temporary cementation


It is not recommended to carry out a temporary fitting of electroformed crowns or bridges. The greatest danger lies in the breaking away of the ceramic facing when the temporary crown is removed. Fixed electroformed ceramic restorations must be cemented permanently.

### ■ Permanent cementation

The established zinc phosphate cements are best suited for permanent cementation. However, modern glass ionomer cements can also be used as well as light curing composite cements (in these cases it is essential to ensure that the working area is relatively or absolutely dry by using cotton wool rolls or similar). The use of zinc phosphate cements and glass ionomer cements is preferable to that of composites when fitting electroformed ceramic restorations. Be sure to use cements with as fine a grain as possible in order to guarantee that electroformed crown keeps its accurately fitting shape even after cementing. It is important to follow the instructions for use issued by the respective cement manufacturers.

It is also imperative that there is an adequate gap for the cement.

Before the restoration is finally cemented into place, the inner surfaces of the crowns must be degreased with alcohol. The tooth dies should be cleaned with a suitable paste and then sprayed carefully with water. Shortly before the restoration is fitted it should be dried again with a blower. **Under no circumstances should the crowns be completely filled with cement.** Because of the excellent fit of the electroformed framework it is sufficient to apply a thin coat of cement to the whole of the inner surface of the restoration with a brush. Then the restoration can be pushed onto the die. In the case of very parallel abutment tooth pre-

parations it may be necessary to assist the cement by using a pumping action in order to push excess cement towards the apical edge of the restoration. When the restoration is seated in its final position, it is advisable to have the patient bite on cotton wool rolls or soft wooden blocks. The patient must hold this position for the duration of the initial curing time (5–10 min).

### ■ Removing excess cement

Excess cement can be removed after about 10–15 min. If zinc phosphate and glass ionomer cements are used, then this is unproblematic for both subgingival and supragingival crown margins.

## ■ Various uses and fitting of AGC® SupraCaps®

| Use                                                                         |                                                                                        | of<br>individually crafted AGC® SupraCaps®                                                                                                                                                                                         |                         |                                                                                                                                                                                              |
|-----------------------------------------------------------------------------|----------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| as                                                                          | using                                                                                  | on modified milled posts<br>or individually crafted<br>implant posts                                                                                                                                                               | on<br>prepared<br>teeth | on individually<br>crafted<br>primary teles-<br>cope crowns                                                                                                                                  |
| AGC®<br>ceramic<br>crowns                                                   | <i>temporary fitting</i>                                                               | !!! Not recommended !!!<br>(neither with nor without cement)                                                                                                                                                                       |                         |                                                                                                                                                                                              |
|                                                                             | <i>permanent cementation</i>                                                           | ZnO phosphate cement or glass ionomer cement<br>(temporary cement may be used on restoration with a limited life span)                                                                                                             |                         |                                                                                                                                                                                              |
|                                                                             | <i>permanent adhesive fitting</i>                                                      | Composite cements Metal and prepared tooth must be treated with surface conditioners before fitting                                                                                                                                |                         |                                                                                                                                                                                              |
|                                                                             | <i>screw fitting</i>                                                                   | Horizontal screws are recommended                                                                                                                                                                                                  |                         |                                                                                                                                                                                              |
| AGC®<br>composite<br>crowns                                                 | <i>temporary fitting</i>                                                               | !!! Not recommended <i>without cement</i> !!!<br>possible <i>with cement</i> providing object is not heated above 200 °C                                                                                                           |                         |                                                                                                                                                                                              |
|                                                                             | <i>permanent cementation</i>                                                           | ZnO phosphate cement or glass ionomer cement<br>(temporary cement may be used on restoration with a limited life span)                                                                                                             |                         |                                                                                                                                                                                              |
|                                                                             | <i>permanent adhesive fitting</i>                                                      | Composite cements Metal and prepared tooth must be treated with surface conditioners before fitting                                                                                                                                |                         |                                                                                                                                                                                              |
|                                                                             | <i>screw fitting</i>                                                                   | Horizontal screws are recommended                                                                                                                                                                                                  |                         |                                                                                                                                                                                              |
| Anchor<br>crowns for<br>fixed<br>bridges                                    | <i>temporary fitting</i>                                                               | !!! Not recommended !!!<br>(neither with nor without cement)                                                                                                                                                                       |                         |                                                                                                                                                                                              |
|                                                                             | <i>permanent cementation or adhesive fitting</i>                                       | ZnO phosphate cement or glass ionomer cement.<br>(temporary cement may be used on restoration with a limited life span)<br>and composite cements Metal and prepared tooth must be treated with surface conditioners before fitting |                         |                                                                                                                                                                                              |
| Anchor<br>crowns for<br>semi-<br>removable<br>bridges                       | <i>screw or latch type retainer</i>                                                    | Horizontal screws are recommended                                                                                                                                                                                                  |                         | Possible with horizontal screws                                                                                                                                                              |
| Secondary<br>telescope<br>crowns<br>for removable<br>bridges or<br>dentures | <i>telescope meso-structures (primary: primary crown or post secondary: SupraCap®)</i> | Parallel sides<br>(Max. angle of 2°)<br>determine desired retainer force;<br>glue into the supra structure with chemically hardening cement<br>eg. AGC® Cem<br>(either on model or in mouth)                                       |                         | Parallel sides<br>(Max. angle of 2°)<br>determine desired retainer force;<br>glue into the supra structure with chemically hardening cement<br>eg. AGC® Cem<br>(either on model or in mouth) |


## Technical procedures

### ■ Crowns and anchor crowns

Study the following instructions carefully. Best results are obtained if all times indicated are observed and the recommended materials are used.

#### ■ Model making

Pour a master model in the usual way from the dental impression. Before duplicating check the die for undercuts and cavities and then block them out using the AGC® blocking out wax. A spacer lacquer must be applied to the prepared dies in order to provide sufficient space for the fixation cement!

#### ■ Note:

Do not block out the plaster die with light curing acrylic and do not treat with die hardener, cyanacrylic glue, grease, oil or similar materials. This can lead to reactions with the duplicating material and result in an inaccurate fit.

#### ■ Duplicating

Shake the AGC® hardener well before use.

Dose: Mix the required amount of AGC® Dubli-Gum and hardener at a ratio of 9 : 1.

Working time: 6 min at 23 °C.  
Single duplication is advisable whenever possible.

Hardness: Shore A 24 to 26.

#### ■ Mixing quantities for the AGC® duplicating system

Duplicating ring size 1:

27 g silicone: 3 g AGC® hardener

Duplicating ring size 2:

54 g silicone: 6 g AGC® hardener

Duplicating ring size 3:


108 g silicone: 12 g AGC® hardener

Duplicating ring size 4:

162 g silicone: 18 g AGC® hardener

Duplicating ring size 4 is intended for denture bases and models.

To facilitate deflasking of the mould, line the metal ring first with foil. (After setting, you can work without the metal ring).

Now briefly stir the duplicating paste and hardener with a spatula and mix for 40 seconds under vacuum. Pour the material in a thin stream from a height of approx. 30 - 40 cm next to the dry die in the flask. Working time approx. 6 min. After 30 min you can remove the die from the mould with the aid of compressed air.

#### ■ Note:

Used AGC® Dubli-Gum is environmentally safe and can be disposed of along with domestic waste.

#### ■ Duplicating the die

Spray the AGC® neutralizer into the mould, allow it to take effect briefly, and then carefully blow dry.

The duplicating mould must be completely dry.

Mix AGC® super hard plaster and distilled water in the following proportions: (100 g AGC® super hard plaster, 20 ml distilled water.)

Mix the AGC® super hard plaster under vacuum for 30 seconds and pour into the mould on a vibrator, ensuring that it is free of bubbles.

Working time approx. 8 min (at room temperature).

After approx. 25 min (from the time of casting) the duplicate die can be removed from the mould. Allow to dry for 2 hours at room temperature.

If you are working with a microwave, proceed as described in the section entitled "Microwave drying".

Grind back the base of the die to make it as small as possible. Round off any sharp edges on the base.

#### Note:

At this stage the duplicate die may no longer be corrected with wax, acrylic or any other material. Any manipulation of the duplicate die can lead to poor deposition and result in an inaccurate fit later!

### ■ Microwave drying\*

The plaster duplicates can be dried in the microwave oven depending on their size, for 1 to 2.5 min.

| AGC® super hard plaster | Anterior die            | Posterior die   |
|-------------------------|-------------------------|-----------------|
|                         | Approx. 1.5 g           | Approx. 2 g     |
|                         | Approx. 2 min – 2.5 min | 1.5 min – 2 min |

You can use any commercially available microwave oven with a rotating glass plate which can be set to a power level of approx. 200 W.

Arrange the parts (not more than 12 dies or 2 - 4 large parts) symmetrically around the center of the rotating glass plate to ensure homogenous microwave drying.

Then remove the dies carefully from the oven, allow them to cool for 1-2 minutes depending on the temperature, and then finish them as usual.

\*"Microwave drying of dental plaster models", Dr. J. Laubersheimer, Quintessence 4/1998.

### ■ Fitting the AGC® contact rod

The procedure for fitting the contact rod is described in the operating instructions accompanying the various models of AGC® electroforming unit.

### ■ Applying the AGC® conductive silver lacquer

Shake the AGC® conductive silver lacquer well before use! You should be able to hear the mixing ball when shaking. The plaster die must be dry and dust-free!

Metal surfaces to be directly electroformed must be grease-free.

Apply the AGC® conductive silver lacquer evenly up to the preparation margin. Dilute the AGC® conductive silver lacquer only as necessary, and take care not to overdilute. Do not apply more than one coat. Correct any defective areas only by touching up the spots themselves.

When using the AGC® Pen, please follow the accompanying instructions for use.


### ■ Important:

In order to make the connection between the conductive silver lacquer on the surface to be electroformed and the AGC® contact rod, apply a coat of AGC® conductive silver lacquer approx. 1 cm wide all around the contact rod and connect the rod to the surface with a thin lug of lacquer.

When using pre-shrunk contact rods, ensure that the shrink fit tubing begins approx. 1 mm away from the die and that the metal rod showing and the start of the isolated part are both entirely covered with conductive silver lacquer.

### ■ Note:

The AGC® conductive silver lacquer coating on plaster dies must be allowed to dry for **15 minutes** at room temperature in order to avoid the risk of residue from the organic binding agent interacting with the AGC® gold electrolyte. The drying time on metal surfaces is at least **30 minutes**.

### ■ Fitting the AGC® shrink fit tubing

The procedure for fitting the shrink fit tubing is described in the operating instructions accompanying the various models of AGC® electroforming unit.

In the case of the AGC® Speed unit no shrink fit tubing is needed!

### ■ Electroforming the parts

The prepared parts can now be electroformed. Refer to the operating instructions accompanying the various models of AGC® electroforming unit.

### ■ Removing the contact rods

For the AGC® Micro, Micro Plus and 5 Process: First, cut off the copper rod close to the plaster die. Then the plaster can be dissolved by placing in AGC® plaster remover in an ultrasonic unit.

For AGC® Speed: Twist the mounting head to remove the contact rod. Place the electroformed die including the contact rod rest into AGC® gypsum remover to dissolve the plaster. Then remove the gold lug by carefully cutting through with a separating disk.

For units with magnetic contact twist the plaster die slightly to remove. The contact rod is reusable.

### ■ Removing the layer of conductive silver lacquer:

If possible, the conductive silver lacquer should be boiled out using unused, diluted nitric acid (30 %) in an extractor. As an alternative, the conductive silver lacquer can be removed by placing the electroformed parts in an acid resistant plastic beaker in an ultrasonic at approx. 60 °C (for approx. 30 min). If acid is not used, the plaster can be removed by sandblasting with clean aluminium oxide (**approx. 110 µm**) at a **maximum pressure of 2 bar**. Examine the surfaces meticulously to ensure that all the conductive silver lacquer has been removed, this is to avoid any discoloration from occurring later on.

### ■ Note:

Conductive silver residues can cause distortions at the margin of the crown and/or discoloration of the ceramic material, and in exceptional circumstances can lead to a discoloration of the gingiva. In the case of secondary parts of telescopic crowns, the conductive silver should **not** be blasted out, since the internal telescopic surface would become roughened, which could result in a poor fit.

In this case the layer of conductive silver lacquer must be boiled out in nitric acid.

### ■ Margin correction

Silicone polishers are recommended for grinding down areas which have been electroformed beyond the margin. Do not use burs.



### ■ Ceramic veneering

Adhesion of the ceramic material to the electroformed substructure is not achieved by means of bonding oxides, but by physical forces and mechanical retention. Adequate adhesion is therefore dependant on the mechanical retention of the electroformed surface.

This retention can be achieved by sand-blasting the electroformed surface with clean aluminium oxide (**approx. 110 µm**) at a **maximum** pressure of **2 bar** (see indicator on machine). Since every sand-blasting unit achieves different results depending on the size of the nozzle and pressure, it is advised to try out on a test coping first. (The surface should look like the coping pictured above.)

Afterwards the AGC® coping is steam cleaned and degreased.

**AGC® gold bonder may be used additionally. (Firing temperature 920 °C).**  
**Please follow the instructions for use enclosed with the product.**

Any commercially available metal ceramic can be used (as long as the firing temperature does not exceed 950 °C).

Please note that the ceramic thickness should be no less than 1 mm.

#### ■ Firing temperature

In the case of the following processes, the firing temperatures should not exceed 950 °C (at a climbing rate of 55 °C/min). Please check whether firing times indicated by the manufacturer can be extended if necessary.

#### ■ Caution:

After the first firing do not use artery forceps or diamond tipped reverse action clamps to hold the electroformed coping, due to danger of distortion.

#### ■ Note:

A temporary fitting of the electroformed crown is contra-indicated. (see "dental procedures").


#### ■ Resin veneering

In order to create a resin veneered electroformed crown, the metal surface must first be treated with a special surface conditioner. Avoid heating to over 200 °C.

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Then carry out the resin build-up according to the manufacturers instructions.

## ■ AGC® electroformed bridge – sintering procedure

The sintered electroformed bridge is made using fine gold in a paste form (Keradec fine gold paste) which creates a secure and strong bond between the AGC® fine gold coping and the separately cast bridge substructure. This bonding method is ideal in terms of its material properties and bio-compatibility, since here the pure gold of the AGC® coping is sintered with fine gold to bridge substructures made of high gold content bonding alloys (usually containing over 85 % gold), eg. AGC® Speziallegierung.

### ■ Preparation

The anchorage of the bridge unit is achieved by a ring shaped structure which surrounds the occlusal cusp tips or incisal edges and palatal surfaces of the abutment teeth.


To ensure adequate stability of the bridge unit, the preparation of the abutment teeth is of great significance. This should include a reduction of 1.5 mm tooth substance with a strong circumferential chamfer or shoulder with rounded inner edges and a box preparation in the approximal region towards the bridge structure as well as the opposing side, which can absorb the horizontal and vertical masticatory forces. Furthermore an anatomical and functional form is important in order to support the ceramic statically. This is done best by using an O-shaped framework design.

To do this fashion the acrylic material from the approximal box on the one side, working outwards in a buccal and oral direction, extending beyond the occlusal edge and into the box of the opposite approximal surface at the end in an O-ring shape. This design enables the possibly lacking anatomical form to be compensated for with small cusp tips, thereby creating an anatomical, ceramic supporting occlusal surface once again.

The metal alloy occlusal rests must not be less than 0.3 mm thick; the material thickness of the ceramic should be between 1 and 1.2 mm.

### ■ Production of the anchor crowns on the abutment teeth

First electroform the copings to a thickness of approx. 0.2 mm (fig. 1). Now transfer the electroformed copings (after finishing and having checked the fit on the master model) to the unsawn working model. Smooth the surface of the gold copings in the area to be worked upon using a rubber polisher (fig. 2). Now, with the use of modelling acrylic, an absolutely stress-free bridge substructure can be produced (no stress to the bridge substructure due to movable sawn dies!).

### ■ Production of the substructure

First apply the acrylic with a small brush using the powder-liquid technique in a ring on the occlusal surface (fig. 3). In order to prevent fracture when lifting off and to make sure that the metal can later flow freely, a material thickness of 0.5 mm



is recommended for the acrylic anchorage, which can then be reduced to 0.3 mm after casting. The individual anchorages can now be connected stress free with a small acrylic connecting bar. The pontic can then be completed in wax (fig. 4). To guarantee an adequate supply of molten metal to the pontic, a pear shaped sprue with a volume comparable to that of the pontic should be attached. After casting the substructure is fitted and thinned out. The cast bridge structure must sit without any gaps on the electroformed copings.

The actual connection of the electroformed copings with the skeletal O-shaped bridge structure is made as follows:

First de-gas the cast substructure according to the alloy manufacturer's instructions.

Then clean and roughen the surface by sandblasting with 110 µm aluminium oxide at a pressure of 2 bar. The electroformed copings are also sandblasted (110 µm max. 2 bar pressure), in order to guarantee optimum sintering with Keradec (fine gold in the form of a paste). The actual sintering process is carried out on a previously made refractory model, which was produced by casting the impression once more with a refractory die material. Now check the fit of the electroformed copings and then the substructure. The entire construction must have a marginally accurate and stress free fit (fig. 5).

The next working step is to remove the bridge structure once more and apply a generous coating of Keradec to the inner surface of the crown anchorage (fig. 6). A thin layer Keradec should also be applied to the area of fixation on the electroformed crown. Now, using slight pressure, place the substructure once more on the electroformed copings on the refractory model and remove excess Keradec with an instrument.


The subsequent drying process is of great significance, since all organic constituents of Keradec must be removed, in order to guarantee that connections and build-up are bubble free.

Dry the object in a ceramic furnace with the firing chamber open at a temperature of approx. 180 °C. The closing time is approx. 10 minutes. The heating up time is 55 °C/min to the end temperature of 940 °C. The end temperature must be held for 4 minutes without vacuum in order to achieve an optimum degree of sintering (fig. 7).

In order to achieve an homogenous, golden-yellow shade of electroformed copings and cast substructure, the entire framework is coated with the AGC® gold bonder (fine gold paste with ceramic constituents) (fig. 8). The gold bonder should not only be painted on, but massaged in well using a short haired brush.

Follow the firing instructions for the AGC® gold bonder.

After having checked the fit, the bridge can be built-up with ceramic. Check that the ceramic you are about to use matches the CTE of the cast gold alloy substructure.

## ■ AGC® SupraCaps® on implant posts


In the field of suprastructures, the individually made SupraCaps® are simple and rational to produce with the aid of the electroforming technique. AGC® SupraCaps® can be made directly and with greatest precision no matter which implant system is used.

The AGC® SupraCaps® on implants are suitable for restorations which are later cemented, screwed or use the telescope technique. The greatest advantage in the case of larger restorations lies in the stress free, passive fit (for fitting instructions see table on page 9).

Select a suitable construction post according to the tooth to be replaced and the implant system used (fig. 1). In the case of preparable gold-titanium posts, the cervical contour of the preparation margin is adapted to the anatomical contour of the gingiva. In the case of cementable restorations, the preparation margin should be situated at a maximum of 1 mm subgingivally.


The industrially prefabricated post shown above (Wieland, BioPost) needs only slight anatomical adaption (figs. 2, 3). It is best to grind the post on a laboratory implant. Cross-cut tungsten carbide burs are recommended. The surface should be worked upon first with coarse and then with fine rubber polishers. The reduction should follow the collar-shaped gingival contour. The cross-section of the post in the area of the soft tissue is adapted to

the cross-section of the root of the tooth to be replaced. Due to the overdimensional size of the industrially manufactured post, it is possible to grind and achieve the optimum gingival contour. The retentive surface is likewise designed according to the anatomy of the tooth. In the case of anteriors, care must be taken in particular to reduce the palatal surface sufficiently. Likewise, the anatomical curvature of the labial surface must be taken into account when preparing the posts, in order to allow for an even layer thickness of ceramic later.

The AGC® SupraCap® is deposited in a direct procedure. The finished prepared post is screwed onto a laboratory implant and the base then invested in pattern resin. Here a die shaped acrylic block should be created, in which the laboratory implant is completely embedded (fig. 4). Ensure that the acrylic is applied directly up to the prepared margin of the die. The screw canal must also be filled with pattern resin. The preparation for deposition and adaption of the contact rod is the same as for the AGC® telescope crown technique.

#### ■ Important:

The hole for the contact rod must be at an adequate distance from the post (fig. 5). The conductive silver lacquer must be applied evenly and cover the surface completely, in order to prevent direct metal connection between the AGC® coping and the post. Because of the metal surface, the conductive silver lacquer must be allowed to dry for at least 30 minutes. In this case the conductive silver lacquer is used as a spacer as well as for separating purposes.

After deposition (fig. 6), the excess electroformed gold is reduced at the preparation margin using an abrasive rubber polisher (fig. 7). The SupraCap® can now be simply removed from the post with the aid of the remaining laboratory implant.

The conductive silver residue in the inner surfaces of the cap must be removed in nitric acid.

#### ■ Important:


Do not remove the conductive silver lacquer by sandblasting. This would damage the finished, marginally accurate gold margin. Remove the pattern resin from the laboratory implant using warmth.

When sandblasting, bonding and building-up the AGC® SupraCap® with ceramic, proceed as for the AGC® single crown.

The final marginal corrections are made on the post. Then, the delicate electroformed margin and the outer surface of the implant post are both polished to a high lustre. These working steps are also achieved with the aid of a previously mounted laboratory implant. This protects the basal area of the post from damage and makes handling easier.

#### ■ Important:

Due to the minimal cement gap, an extremely fine grained and runny cement must be used for fixation! (fig. 8, 9)

## ■ AGC® telescope crown technique

### ■ Producing the primary parts

The primary copings are conventionally milled in wax. We recommend parallel milling (telescope crowns) with a circumferential shoulder which stabilizes the later tertiary construction as well as the ceramic build-up. A small gold margin should be left at the shoulder, showing a slightly overcontoured anatomical form (fig. 1).

### ■ Important:

After the secondary parts have been electroformed, the excess gold is reduced at the margin by simultaneously grinding the shoulder. This also happens when contouring the tertiary structure after glueing. For this reason, the shoulder should be designed somewhat broader, in order to have sufficient material later for the anatomical form.

### ■ Exception:


When designing a cover denture, do not include a circumferential chamfer in order to allow for resiliency.

The primary parts are made according to conventional casting procedures. A copper free alloy with a high gold content must be used. In the case of alloys contain-

ing copper, there is a risk of discoloration (darkening) of the surface due to copper oxide formation.

For the primary crowns we recommend the copper free AGC® Speziallegierung (Wieland Dental+Technik).

After removing the sprue and checking the fit on the master model (pay attention to air bubbles and an exact fit) a milling die is made. Then place all primary parts on the milling model.

### ■ Important:

First roughly pre-mill the surface using a coarse, cross-cut, parallel bur (fig. 2). Then smooth the parallel surface more and more using fine burs and milling oil. Optimal functioning of the telescope technique through adhesion depends almost entirely on the surface smoothness of the primary crown (fig. 3).

For final finishing we recommend using a parallel bur, wrapped a light layer of cotton wool and diamond polishing paste. It is imperative that the primary crowns are milled and polished in a milling machine. If the crowns are polished by hand it can easily create an uneven surface which means that the gold secondary crown deposited directly onto the primary crown, is impossible to remove.

The previously made circumferential shoulder makes it easier to separate and remove the crowns after electroforming, and gives the electroformed secondary crown sufficient stability in the marginal area. Optimum retention is achieved when the parallel sides are as long as possible. In the case of anterior teeth a compromise should be found between palatal sloping and the circumferential parallel surfaces.

#### ■ Preparation for the direct electroforming of the secondary parts

Fill the polished primary crowns with pattern resin by placing the telescope crown upside down onto the duplicating system's red fixation wax. The circumferential margin must be completely covered with resin up to the edge of the shoulder so that no gold can be deposited there during electroforming (fig. 4). The circumferential margin should now be covered in acrylic approx. 1 mm thick. The resin must not taper to a thin end,

since otherwise it could peel back in the hot gold bath. After setting the margin is checked and adjusted if necessary.

By including a stainless steel screw into the pattern resin, the primary crown can be held better later therefore easing removal of the deposited secondary crown. The screw must be removed before electroforming and can be replaced afterwards.

The contact rod is attached 1-2 mm underneath the crown margin.

Before applying a thin, even layer of conductive silver lacquer, the metal surface must be de-greased with alcohol (AGC® de-greaser, Wieland Dental+Technik). Shake the AGC® conductive silver lacquer well before use! (You should be able to hear the mixing ball when shaking). Apply a thin but still covering layer to the surfaces to be electroformed. Make sure the layer is very thin. The AGC® conductive silver lacquer acts as a spacer for the saliva later. The thinner the conductive silver lacquer, the narrower the gap and the better the adhesive power. There must however, be absolutely no holes in the lacquer layer.

Please do not forget the connection between the conductive silver lacquer and the contact rod (fig. 5).

The AGC® conductive silver lacquer must be allowed to dry for at least 30 minutes.

## ■ Electroforming the secondary crowns

Please see the appropriate instructions for use for the various AGC® electroforming units by Wieland Dental + Technik.

## ■ Important:

Please observe the following when electroforming the secondary crowns in the AGC® Speed machine!

To achieve a deposition thickness of 200 µm, enter 300 µm on the display. Then, after choosing the appropriate power level, enter one power level below.

This is not relevant for a thickness of 300 µm.

## ■ Further processing of the electroformed secondary crowns

After the deposition time has elapsed, the crowns can be removed from the AGC® unit. First remove the contact strip. Then pull out the contact rod (fig. 6).

In order to hold the die more firmly, replace the stainless steel screw into the pattern resin.

If the cervical margin was well isolated, the primary and secondary crowns can already be easily removed. Should this not be the case, the electroformed margin must be reduced with a silicone polisher down to the circumferential shoulder of the primary crown (fig. 7).

## ■ Important:

If the distance from the crown margin to the milled margin is very small, please reduce very carefully in order not to damage the crown margin. When reducing the AGC® margin with an abrasive rubber polisher, the angle should once again be anatomically overcontoured (fig. 8). Should the secondary crown be difficult to remove, little, individually made, pattern resin handles on the outside can aid in pulling off (fig. 9) or a riving hammer held against the screw may help. The vibrations help to loosen the secondary crown. A piece of rubber dam can give a better grip when pulling the crowns apart.

After the primary crowns have been separated from the secondary crowns, the conductive silver lacquer on the inner surfaces must be removed. This is achieved by boiling the secondary crown in nitric acid (30 %). As an alternative, the

conductive silver lacquer can be removed by placing the electroformed parts in an acid resistant plastic beaker with fresh nitric acid in an ultrasonic at approx. 60 °C for approx. 30 min. Then discard the nitric acid and use a fresh amount the next time.

#### ■ Caution:

Always boil under a closed extraction unit! Open the acid container under the extraction unit. Follow the safety instructions!!

#### ■ Important:


The conductive silver must be completely boiled out, otherwise this could cause black discolorations on the inner surfaces of the telescopic crown due to oxidation of the conductive silver lacquer. The conductive silver must never be sandblasted out in the AGC® telescopic crown technique, so that the inner surfaces are not roughened, which would prevent the gli-

ding effect (fig. 10). Any residue after boiling out is removed with a steam cleaner or an ultrasonic.

#### ■ Production of the Tertiary supra-structure for affixing

Secondary parts manufactured by electroforming can be affixed to tertiary model casting structures, removable bridges, implant suprastructures and cover denture prostheses.

Since the secondary crown is later inserted and glued into the tertiary construction, space must be left for the adhesive. The layer thickness necessary to achieve complete hardening of the adhesive is approx. 0.1 mm. In the case of tertiary structures, space is created by means of spacer lacquer before duplicating. For removable bridges in the direct procedure, an adaptafoil (fig. 11) with a corresponding spacer is used, which is of course removed before casting.

## The tertiary structure

It is generally the case that the electroformed secondary part is partially or completely covered by the tertiary construction depending on the indication. This must always be decided according to the static aspects of each individual case. The metal substructure must be constructed in such a way that the secondary parts are at least supported in a claw shaped manner mesially and distally with an occlusal support. If there is enough space, the secondary parts can also be circumferentially supported and provided with a protective plate at the back.

It is necessary to make sure of sufficient adhesive areas, eg. protective back plate and metal occlusal surfaces.

## Glueing

The connection between electroformed secondary parts and the suprastructure is achieved by means of the adhesive technique. This technique has clear advantages compared with other techniques for joining parts together:

§Stress free connections

§Simple handling

§Heat free technique – therefore no structural change of the parts to be joined.


In order to obtain an optimum bond between the parts to be joined and the adhesive, the surfaces must first be conditioned. Use 110 µm aluminium oxide with a max. pressure of 2 bar and metal primer II (by GC) or AlloyPrimer (Kuraray). Please read the working instructions.

Never use conditioning systems which need heat to activate (loss of hardness).

When both parts to be joined have been conditioned accordingly, they are glued with a self activating compomer cement. We recommend the specially developed AGC® Cem (Wieland Dental + Technik), which due to esthetical reasons has the color of honey (fig. 12). The primary

crowns can be protected with vaseline.

Place the secondary crowns carefully onto the primary crowns. Do not touch the conditioned parts with your fingers. Make sure the fit is accurate. Apply the two-component adhesive to the electroformed secondary crowns as well as to the tertiary structure. Affix stress free on the master model. After hardening, take any excess glue off and remove the complete restoration from the master model. Now the final anatomical form around the circumferential adhesive gap can be trimmed and polished.

## ■ AGC® Bar sleeves

### ■ Preparation

Construct the individually milled bar in the usual manner.

Bar sleeves should be made on bar attachments which have either been milled parallel or at a maximum of 2°!

### ■ Note:

The bar attachment must be thoroughly cleaned, free of polishing material and other residues (ultrasound cleaning). There must also be no undercuts.


### ■ Preparing the bar attachment for electroforming

Place the bar onto the master model (fig. 1). With implants in the end position screw the bar securely on to the master model, block out the interdental spaces and beneath the bar. Ensure there are no undercuts.

Duplicate the blocked out model with the bar in position. For small bar constructions the silicone matrix duplicating

method may be used. Remove the model. Remove the bar and clean the wax residues off (using boiling water or a steam cleaner).

Fill abutment crowns with pattern resin the same as in the telescope crown technique and insert a stainless steel screw into the resin (fig. 2). Remove the stainless steel screws before electroforming. After electroforming the screws are replaced to help give a firm grip when removing the bar sleeve. If the bar contains implants the laboratory analog is screwed in place and the screw canal is sealed with pattern resin.

The bar with the laboratory analogs and/or screws is replaced in to the duplicated form and cast with Alpa Pur (Alpina) model acrylic or with AGC® Special Plaster.

The model implants must be completely covered with pattern resin to avoid gold being deposited onto them during electroforming. The stainless steel screws are removed before electroforming.

### ■ Surface area and gold calculation

Bar attachments may be electroformed in the AGC® Micro with a surface area of up to 1500 mm<sup>2</sup> and a thickness of 0.2 mm.

With a thickness of 0.3 mm the surface area must not exceed 1120 mm<sup>2</sup>. It is advisable to electroform a thickness of 0.3 mm.

The AGC® Micro Plus has no restrictions.

First the total surface area must be calculated in mm<sup>2</sup>. Choose the power levels of the abutment teeth or the implants from the comparison chart. Then calculate the surface area of the structure to be electroformed in mm<sup>2</sup>.

■ The formula for calculating the bar surface is as follows:

Surface

$$F = (\text{width} \times \text{length}) + 2 \times (\text{height} \times \text{length})$$

The power levels (Micro / Micro Plus) are equivalent to the following values:

Level 1 ≈ 50 mm<sup>2</sup>

Level 3 ≈ 100 mm<sup>2</sup>

Level 5 ≈ 160 mm<sup>2</sup>

Level 7 ≈ 250 mm<sup>2</sup>

Level 2 ≈ 70 mm<sup>2</sup>

Level 4 ≈ 120 mm<sup>2</sup>

Level 6 ≈ 200 mm<sup>2</sup>

Example 1: (fig. 3)

2 abutment teeth (level 4 and level 5) plus bar sleeve (length 25 mm, height 4 mm, width 3 mm).

The total surface area is calculated as follows:

2 abutment teeth =

$$120 \text{ mm}^2 + 160 \text{ mm}^2 = 280 \text{ mm}^2$$

bar sleeve =

$$(3 \text{ mm} \times 25 \text{ mm}) + 2 \times (4 \text{ mm} \times 25 \text{ mm}) = 275 \text{ mm}^2$$

The total surface area is

$$280 \text{ mm}^2 + 275 \text{ mm}^2 = 555 \text{ mm}^2$$

Power level 7 is needed twice and power level 1 is needed once. At a thickness of 0.3 mm, 258 ml gold electrolyte will be needed. 2 contact rods must be attached to the object. Each rod will need 4 ml of brightening additive.

We recommend using the AGC® Spider for setting the contacts.

Example 2: (fig. 4)

Bar attachment on 4 implant abutments, length 84 mm, width 2 mm, height 6 mm

Surface F =

$$(\text{width} \times \text{length}) + 2 \times (\text{height} \times \text{length})$$

$$\text{Result: } 1176 \text{ mm}^2$$

Power level 7 is needed 4 times and power level 5 is needed once.


At a thickness of 0.3 mm, 540 ml gold electrolyte will be needed. As the AGC® Micro has a gold electrolyte limit of 523 ml, the power level 5 must be changed to power level 4 in order to reduce the amount of liquid needed to 523 ml, the maximum fill level. 3 contact rods must be attached. **Each** rod will need 4 ml of brightening additive.

We recommend using the AGC® Spider for setting the contacts.

With a chosen thickness of 0.3 mm the bar sleeve result will be slightly thinner.

#### ■ Important:

If the gold electrolyte volume in the AGC® Micro unit exceeds 450 ml, the magnetic agitator rod for bar attachments must be used. (AGC® bar magnetic agitator rod).

5

6

#### ■ Connecting the parts

The contacts must be glued as shown in figs. 5 and 6.

Contacts must not be applied to the ends of the bar.

Before applying the AGC® conductive silver lacquer to the area to be electroformed, all loose pieces of plaster or pattern resin must be removed using the compressed air and then the surface must be thoroughly degreased, with AGC® degreaser.

AGC® conductive silver lacquer is applied in an even thin layer up to the end of the milled surface / margin and allowed to dry for at least 60 minutes. Then apply the conductive silver lacquer in a wide connection between the bar and the contact rods (3 mm).

To create a clean, smooth joint between bar sleeve and acrylic, electroform the bar sleeve 1-2 mm onto the alveolar ridge. This means the conductive silver lacquer must be applied 1-2 mm over the alveolar ridge.

### ■ Caution:

This must be taken into consideration when calculating the surface area.

### ■ Electroforming

The bar attachment to be electroformed should be placed into the AGC® electroforming machine as follows:

The surface to be electroformed faces the outside of the glass, away from the anode. The bar ends face the base of the glass (figs. 7 + 8).

### ■ Separating and removing the bar sleeve

Apply little pattern resin handles onto the outside of the galvanised bar sleeve to act as removal aids. With help of the removal aids and a riveting hammer the bar sleeve can be easily parted from the bar. The prerequisite is of course a perfectly milled surface.

The model implants or the stainless steel screws act as anchorage elements during separation.

### Tertiary structure

A tertiary structure is absolutely necessary for strength and stability. One retention around the bar sleeve is enough. The AGC® bar sleeve is glued / adhesively bonded into this retention.

### Glueing


We recommend the AGC® Cem for adhesive bonding (fig. 9). The surfaces to be glued should be treated with Metall-primer II first.

Afterwards the piece of work can be finished in the usual manner.

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© WIELAND DENTAL + TECHNIK
Dentalite, S.A.
Chin Min Tai Enterprise Co. Ltd.
Wieland Dental Systems, Inc.